


PulsePoint

ACTIVATE YOUR COMMUNITY

The chain of survival is only as strong as we make it. Create a culture of action in your community.



Richard Price
Founder/President
PulsePoint Foundation



PulsePoint.org | PulsePoint is a 501c3 Non-Profit Foundation



SCA victim in need.





**911 system sends
PulsePoint alert.**




**Signal received by nearby
PulsePoint users.**



**Users rush to help
the victim before
professional help arrives.**





SYSTEM PERFORMANCE

2m 54s

EARLY ARRIVAL ADVANTAGE ¹

5m 07s

COMMUNITY NOTIFICATION
OPPORTUNITY ²

1m 52s

COMMUNITY RESPONSE TIME

INCIDENTS

40,924

INCIDENTS LAST 24 HOURS

864

SCA INCIDENTS LAST 24
HOURS ³

1,275,743

SCA INCIDENTS TO DATE

PulsePoint Report

Utilizing CARES 2018 - 2020 Data

2018
 Total # of PulsePoint cases in CARES = 18,361
 Total # of PulsePoint Agencies participating: 168/371 (45.3%)

2019
 Total # of PulsePoint cases in CARES = 28,227
 Total # of PulsePoint Agencies participating: 350/510 (68.6%)

2020
 Total # of PulsePoint cases in CARES = 42,242
 Total # of PulsePoint Agencies participating: 407/574 (70.9%)

	Non-PulsePoint* 2018	PulsePoint 2018	National 2018
Non-Traumatic Etiology Survival Rates			
Overall	10.0% (63503)	11.6% (18361)	10.4% (81864)
Bystander Wit'd	15.2% (23884)	18.2% (6685)	15.8% (30569)
Unwitnessed	4.2% (31643)	5.1% (9406)	4.4% (41049)
Utstein ¹	32.1% (6860)	37.1% (2186)	33.3% (9046)
Utstein Bystander ²	36.1% (3919)	40.7% (1386)	37.3% (5305)
Bystander Intervention Rates³			
CPR	37.8% (47098)	43.6% (13624)	39.2% (60722)
Public AED Use	11.4% (8315)	13.4% (2622)	11.9% (10937)

	Non-PulsePoint* 2019	PulsePoint 2019	National 2019
Non-Traumatic Etiology Survival Rates			
Overall	10.3% (72729)	11.2% (28227)	10.5% (100956)
Bystander Wit'd	15.5% (27761)	17.2% (10583)	16.0% (38344)
Unwitnessed	4.2% (35811)	4.9% (14196)	4.4% (50007)
Utstein ¹	32.3% (8099)	35.3% (3372)	33.2% (11471)
Utstein Bystander ²	36.4% (4799)	39.4% (2149)	37.3% (6948)
Bystander Intervention Rates³			
CPR	39.8% (54017)	44.7% (21110)	41.2% (75127)
Public AED Use	11.3% (9582)	14.3% (3971)	12.2% (13553)

	Non-PulsePoint* 2020	PulsePoint 2020	National 2020
Non-Traumatic Etiology Survival Rates			
Overall	8.6% (85134)	9.6% (42242)	9.0% (127376)
Bystander Wit'd	12.7% (31628)	14.6% (15578)	13.3% (47206)
Unwitnessed	3.9% (42477)	4.6% (21343)	4.1% (63820)
Utstein ¹	27.7% (8147)	32.1% (4094)	29.2% (12241)
Utstein Bystander ²	31.7% (4722)	35.5% (2547)	33.0% (7269)
Bystander Intervention Rates³			
CPR	38.4% (62853)	43.8% (31840)	40.2% (94693)
Public AED Use	8.6% (8865)	9.9% (4343)	9.0% (13208)

An analysis of PulsePoint agencies who were also CARES participants in 2018, 2019, and 2020 showed that PulsePoint communities performed better in all of the above indicated CARES performance measures. In 2020, bystander CPR was 5.4% higher (43.8% vs. 38.4%) and public AED use was 1.3% higher (9.9% vs. 8.6%) in PulsePoint communities compared to non-PulsePoint CARES communities. Overall survival for sudden cardiac arrests that were witnessed and shockable (Utstein) was 4.4% higher (32.1% vs. 27.7%). In 2020, 70.9% of PulsePoint agencies participated in CARES and 33.2% of CARES cases were from PulsePoint agencies (42,242 of 127,376).

-Please Note: Descriptions of, references to, or links to other products, publications, or services does not imply endorsement of any kind unless expressly stated by CARES. The reports are intended for educational, research, self-study, and reference purposes only and are not intended to substitute for the exercise of professional judgment. CARES makes no warrant to the accuracy, completeness, suitability, or result obtained from the content.

*Non-PulsePoint rates exclude PulsePoint agencies

¹ Utstein: Witnessed by bystander and found in shockable rhythm.
² Utstein Bystander: Witnessed by bystander, found in shockable rhythm, and received some bystander intervention (CPR and/or AED application).
³ Bystander CPR rate excludes 911 Responder Witnessed, Nursing Home, and Healthcare Facility arrests. Public AED Use rate excludes 911 Responder Witnessed, Home/Residence, Nursing Home, and Healthcare Facility arrests.



2020 GUIDELINES

NEWS FROM THE AMERICAN HEART ASSOCIATION

Significant New and Updated Recommendations

Systems of Care

1. *“Summoning willing bystanders: Emergency dispatch systems should alert willing bystanders to nearby events that may require CPR or AED use through mobile phone technology (Class 1, LOE B-NR).”*
2. *Notification of lay rescuers via a mobile phone app results in improved bystander response times, higher bystander CPR rates, shorter time to defibrillation, and higher rates of survival to hospital discharge.”*

Released 10/21/2020

10 STEPS

to Improve Cardiac Arrest Survival



1. Establish a cardiac arrest registry
2. Begin Telephone-CPR with ongoing training and QI
3. Begin high-performance EMS CPR (HP-CPR) with ongoing training and QI
4. Begin rapid dispatch
5. Measure professional resuscitation using the defibrillator recording (and voice if possible)
6. Begin an AED program for first responders, including police officers, guards, and other security personnel
7. Use smart technologies to extend CPR and public access defibrillation programs to notify volunteer bystanders who can respond to a nearby arrest to provide early CPR and defibrillation
8. Make CPR and AED training mandatory in schools and the community
9. Work toward accountability – submit annual reports to the community
10. Work toward a culture of excellence

"The non-PulsePoint dispatch associated OHCA patients experienced a similar rate of bystander CPR to the [CARES] published figure (35.5%), while the PulsePoint dispatch-associated group experienced a significantly increased bystander CPR rate of 59.4%."



ELSEVIER

Available online at [ScienceDirect](#)**Resuscitation**journal homepage: www.elsevier.com/locate/resuscitation**Clinical paper****PulsePoint dispatch associated patient characteristics and prehospital outcomes in a mid-sized metropolitan area**

Tanner Smida^{a,b,*}, Jessica Salerno^a, Leonard Weiss^a, Christian Martin-Gill^a, David D. Salcido^a

^a University of Pittsburgh School of Medicine, Pittsburgh, PA, United States

^b West Virginia University, United States

App Version / Responder Types

Public App

- Who
 - **CPR-trained citizens**
- **CPR-needed Alert**
 - **Public locations only**
- **Public safety activity**
 - **Limited to public information that is useful to community members**
 - Vehicle accidents, utility emergencies, fires, etc.



PulsePoint Respond 12+


Building informed communities
[PulsePoint Foundation](#)

#5 in Medical

★★★★★ 4.8 • 131K Ratings

Free

Screenshots [iPhone](#) [iPad](#)



App Version / Responder Types

Community Verified Responders

- Who
 - **Doctors, Nurses, CERT, Security Guards, CPR Instructors**
- CPR-needed Alert
 - **All locations**
- Public safety activity
 - **Limited to public information that is useful to community members**
 - Vehicle accidents, utility emergencies, fires, etc.




App Version / Responder Types

Professional Verified Responders

- Who
 - **Firefighters, EMT's, Paramedics**
- CPR-needed Alert
 - **All locations**
- Public safety activity
 - **Complete incident information including potential HIPAA data**






PulsePoint Public Responder

Target user: *Community members* interested in staying informed of local emergency activity and those trained in CPR and willing to assist if a cardiac arrest occurs in a public place near them.


- Nearby CPR-needed alerts in **public locations**.
- **Filtered list of incidents** displayed.
 - Calls of a dangerous or sensitive nature typically not displayed.
- Medical Emergency incident types **described generically**.
- Full incident **address not shown** on medical incident types.
- Routing (turn-by-turn) **not offered on medical incident types**.
- Informed community alerts (**wildfires, flooding, etc.**).
- View and add AED locations, including colocated items.
- Self-selected participation.



PulsePoint Verified Responder - Pro

Target user: *Public safety employees*, typically firefighters, paramedics, and law enforcement officers assigned Verified Responder Pro status by their organization.


- Nearby CPR-needed alerts in **all locations**.
 - Digital ID badge
- **All incidents** displayed.
- Expanded medical incident types.
- Full incident **address shown** on medical incident types.
- Routing (turn-by-turn) **offered on medical incident types**.
- Destination Street View (Android)
- Additional notification types (**MCI, Mutual Aid**).
- Add and **approve** AED locations, including colocated items.




PulsePoint Verified Responder - Community


Target user: *Invited members of the community with medical or rescue training such as public safety retirees, CERT members, doctors, nurses, and residential security staff or residents with special neighborhood responsibilities.*

- Nearby CPR-needed alerts in **all locations**.
 - Digital ID badge
- Same filtered incident list and restrictions as the public version.
- Add and **approve** AED locations, including colocated items.



VS.





9:41

Back

Location/Business

Diablo Vista Middle

AED Location Description

Inside labeled wall

Restricted Access

AED is not in cart, etc.)

Street Address

4100 Camino Tassa

City

Danville

County

Contra Costa Coun

Zip Code

94506

State

California

9:41

Back



Drag

Apple Maps

9:41

Back



9:41

Back

Indicate additional contents for this AED.

Bleeding Control

Bleeding Control Expiration

February 20, 2022

Epinephrine

Epinephrine Expiration

June 27, 2021

Naloxone

Naloxone Expiration

September 15, 2021



9:41

Back

Manufacturer

Physio-Control

Model

LIFEPAK CR2

Asset ID

Serial Number

48474132

Install Date

November 20, 2019

Battery Expiration

November 23, 2023

Electrode Expiration

July 24, 2024

Pediatric Electrode

Medical Direction



9:41

Cancel

Remove

Issue *

Business is closed

Comment

This AED is shown closed recently and will be removed from the map.




AED is shown with incorrect name

Business is closed



9:41

Pin



Map Settings

Map

AED Filters

Pending

Has Image

Bleeding Kit



9:41

Done Burleigh Manor Middle School



AED Location

Location/Business Name

Burleigh Manor Middle School

AED Location Description

Mounted On Wall Outside Health Room

Address

4200 Centennial Lane
Ellicott City, Maryland 21042
Howard County, United States

Agency



Howard County Fire and Rescue

Colocated Items

Bleeding Control

Site Coordinator



Oct 04 2021, 23:15:57

5303 Heron Creek Boulevard

North Port, 34287

(315) 735-6119



Close call

What3Words ///reverted.twisted.cafes

Show all

Lat/Long 27.0598069, -82.2212612

DMS 27° 3' 35" N, 82° 13' 17" W

SMS chat

Location timeline

Additional data

Location update 23:15:57


Approximate address:	Class:	Apple
5303 Heron Creek Boulevard,	Uncertainty:	12 m
North Port,	Confidence:	90%
34287	Altitude:	125 m

6 mins Apple

Search Location...

Heron Creek Golf & Country Club

AgencyID	16072
Address	5303 Heron Creek Boulevard
Description	Fitness Center
City	North Port
PlaceType	
AccessCode	
Manufacturer	
Model	
Collocated items:	None
Accessibility:	Publicly Accessible



Zoom to

Heron Creek Golf & Country Club

Heron Creek Golf and Country Club

Morgan Center

Oct 04 2021, 23:15:57

5303 Heron Creek Boulevard

North Port, 34287

(315) 735-6119

Close call

What3Words ///reverted.twisted.cafes

Show all

Lat/Long 27.0598069, -82.2212612

DMS 27° 3' 35" N, 82° 13' 17" W

SMS chat Location timeline Additional data

Location update 23:15:57

Approximate address:	Class:	Apple
5303 Heron Creek Boulevard,	Uncertainty:	12 m
North Port,	Confidence:	90%
34287	Altitude:	125 m

6 mins Apple

Search Location...

Heron Creek Golf & Country Club



Zoom to

Heron Creek Golf & Country Club

Heron Creek Golf and Country Club

Morgan Center

Oct 04 2021, 23:15:57

5303 Heron Creek Boulevard

North Port, 34287

(315) 735-6119

Close call

What3Words ///reverted.twisted.cafes [Show all](#)

Lat/Long 27.0598069, -82.2212612

DMS 27° 3' 35" N, 82° 13' 17" W

SMS chat Location timeline Additional data


Location update 23:15:57

Approximate address:	Class:	Apple
5303 Heron Creek Boulevard, North Port, 34287	Uncertainty:	12 m
	Confidence:	90%
	Altitude:	125 m

6 mins Apple


463

BS



BUILDING INFORMED COMMUNITIES


PulsePoint notifications provide an early heads-up to local threats such as wildland fires, flooding and utility emergencies.

DOWNLOAD THE APP 



LATEST NEWS:



February is Heart Month! Check out our Heart Month Resource Toolkit—full of ideas and resources to engage and encourage your citizens to take care of their hearts, learn CPR, register AEDs, and download PulsePoint. 

LAST 24 HOURS

90

CPR-NEEDED ACTIVATIONS

225

NEARBY RESPONDERS
ALERTED

386

NEARBY AEDS IDENTIFIED

58

AEDS ADDED

2,134

NEW APP USERS

[More Stats...](#)

PulsePoint



Richard Price

richard.price@pulsepoint.org



General Information

info@pulsepoint.org



Website

pulsepoint.org



Twitter

[@pulsepoint](https://twitter.com/pulsepoint)